



LACTOPLUS

The versatile solution for quality silages

The product

Formula: *L. Plantarum* 14D/CSL and *L. Buchneri* 40788 to dominate the fermentation and protect the silage at the opening of the silo

Concentration: synergy of two pre and post opening strains 4×10^{11} UFC/g

Packaging: envelopes of 100g

Benefits

Stability anaerobic and aerobic for counteract the development of molds and yeasts even after opening the silo

Maximize nutritional value and greatly reduces weight loss of silage in the trenches Sanitary Hygienic

Control of the entire supply chain, thanks to the inhibition of putrefactive and nitrogenous substances

Perfect for large trenches size

A low-angle, upward-looking photograph of a cornfield. The green leaves of the corn plants are the dominant feature, filling most of the frame. In the upper right, the bright sun is visible against a clear blue sky, creating a lens flare effect. The perspective makes the plants appear to be reaching towards the center of the image.

SACCO
system
— agrovet

Colture for silage